Protocol for analytical ultracentrifugation sample preparation:

- 1. Provide partial specific volume (v-bar) of the proteins. If v-bar is not available, send complete amino acid sequence of your protein via e-mail. Please type the sequence into the body of the e-mail rather than sending it as an attached file.
- 2. Provide minimum 250 μ L of each sample to be analyzed for equilibrium runs and 1 mL for velocity runs and at least 5 mL of buffer.
- 3. For accurate determination of molecular weight or oligomeric state, provide three different concentrations for each sample. An example of an ideal concentration series is 0.25 mg/mL, 0.50 mg/mL, and 1.00 mg/mL. Ensure that the absorbance at either 230 nm or 280 nm (or, whatever other wavelength) is not less than 0.25 at the lowest concentration, and does not exceed 1.00 at the highest concentration. Concentrations as low as 1 µg/mL are acceptable only if the absorbance at 230 nm is at least 0.30. Please make your own dilutions.
- 4. Try to avoid using DTT, 2-mercaptoethanol, and protease inhibitors. If DTT is absolutely required to maintain sample stability, do not exceed 5 mM concentration. For ATP or ADP try to use less than 0.3 mM.
- 5. The reference buffer should contain all components that are in the sample minus the protein.
- 6. Provide us with the samples, buffers, and filled form.
- 7. Samples will be stored at **4**°**C** and run at **4**°**C**, unless otherwise specified.

Acceptable buffer components:

Urea

1-Propanol 2-Propanol Ammonium chloride Ammonium hydroxide Cadmium chloride Cadmium sulfate Citric acid Cobaltous chloride EDTA disodium Ethanol Formic acid Fructose Guanidine hydrochloride Hydrochloric acid Lactose Lanthanum nitrate Magnesium chloride Magnesium sulfate Mannitol Methanol Oxalic acid Phosphoric acid Potassium bromide Potassium carbonate Potassium dichromate Potassium ferricyanide Potassium iodide Potassium nitrate Potassium PO₄ di-basic Potassium PO₄ mono-basic Procaine hydrochloride Propylene glycol Sodium bicarbonate Sodium bromide Sodium citrate Sodium diatrizoate Sodium hydroxide Sodium molybdate Sodium PO₄ mono-basic Sodium PO₄ tri-basic Sodium thiosulfate Sodium thiocyanate Sucrose Sulfuric acid Trichloroacetic acid Trifluoroethanol

Zinc sulfate

Acetic Acid Ammonium sulfate Calcium chloride Creatinine Ethylene glycol Glucose Inulin Lead nitrate Maltose Nickel sulfate Potassium bicarbonate Potassium chloride Potassium ferrocyanide Potassium oxalate Potassium sulfate Silver nitrate Sodium carbonate Sodium dichromate Sodium nitrate Sodium sulfate Sodium tungstate Tartaric acid

Tartaric acid
Tris (hydroxymethyl)aminomethane

Acetone
Barium chloride
Cesium chloride
Cupric sulfate
Ferric chloride
Glycerol
Lactic acid
Lithium chloride
Manganous sulfate
Nitric acid

Potassium biphthalate
Potassium chromate
Potassium hydroxide
Potassium permanganate
Potassium thiocyanate
Sodium acetate
Sodium chloride
Sodium ferrocyanide
Sodium PO₄ di-basic
Sodium tartrate
Strontium chloride
Tetracaine hydrochloride